


Signals of muscles and peripheral nerves

Michal Novotny


Muscles

- Smooth muscles
 - Slow automatic movements of organs (oesophagus, stomach, bowels)
- Striated muscles
 - Cardiac muscles – automatic movements of heart
 - Skeletal muscles – voluntary movements




Hierarchical structures for muscles

- Muscle can only contract
- Contraction is caused by emission of acetylcholine to the membrane of muscle fiber -> action potential leads to opening of Ca^{2+} canals and leakage of Ca^{2+} ions to sarkoplasma





Sarcomere





Contraction

- ATP arrival causes bending and the myosin head pulls actin and myosin fibers towards each other

Event order during muscle contraction




Muscle control



Motor unit


- Motor unit = spinal motoneuron + axon + all innervated muscle fibers
- The most basic voluntarily activated muscle unit
- Excitation vs. inhibition



Axon of motor neurons extend from the spinal cord to the muscle. There each axon divides into a number of axon terminals that form neuromuscular junctions with muscle fibers scattered throughout the muscle.


Motor unit EMG generation

- Anatomy:
 - Motor units – ventral spinal horn, moto-neuron, axon, muscle fibers
- Physiological model:
 - Summation of impulses arriving from spinal cord activates motor units generating action potentials
 - Very low voltage ($100\mu\text{V}$) but synchronous in many fibers
- Instrumentation:




Motor unit size

Muscle	# neuronal fibers	# muscle fibers	# motor units	# muscle fibers in motor unit
Platysma	1826	27100	1096	~25
Pointing finger muscles	199	40500	119	~340
Lower limbs muscles	965	1120000	579	~1934



Platysma

<https://eluc.kr-olomoucky.cz/verejne/lekce/188>



Hand muscles

<https://fineartamerica.com/featured/9-muscles-of-the-hand-asklepios-medical-atlas.html>




Lower limb muscles

<https://anatomyzone.com/lower-limb/>

Wave types and action potentials


- Needle electrodes
- Usually bi or triphasic
- 3-15 ms in length, 100 - 300 μ V and 6 – 30 Hz
- Shape is given by used electrode



Monophasic
waves

Biphasic
waves


Triphasic
waves




Motor units activation

- Spatial activation
- Temporal activation
- Quantitative characteristics due to the signal complexity


Easy work



Medium work





Hard work




Motor units activation

- Hanemann's rule – motor units are activated sequentially
- Muscle force is dependent on number of activated motor units
- Muscle fibers
 - Type I (slow 110ms, 50%, endurance, activated first)
 - Type II (fast 50ms, force)
 - Type IIa 25%
 - Type IIx 25% (strongest, activated last)






Muscle stimulation

- Indirect
 - Local potential at neuromuscular junction
- Direct
 - Electric signal causes artificial depolarization







EMG and contraction power

- Muscle force is given by motor unit recruitment and frequency of incoming impulses
- Frequency modulation
- Similar force can be produced by recruitment of smaller MU stimulated by higher frequency or larger MU with lower frequency



EMG and contraction power




EMG recording

- Surface electrodes
 - Less invasive
 - Easy application
 - 20-500Hz
 - Generalized signal

Bagnolli - surface



Nihon-Kohden
Surface, Disc




EMG recording

- Needle electrodes
 - Invasive
 - Up to 5000 Hz
 - Localized
 - Expensive and complicated application


Nihon-Kohden
concentric



Nihon-Kohden
needle



Surface electrode placement




Electromyography

- Applications
 - Diagnostic
 - Needle EMG
 - Conduction study
 - Kinesiology
 - Fatigue analysis
 - Prosthetics
 - Others




Native EMG

- Acquisition during muscle relaxation
- Needle
- Comparison of features of motor units (duration, amplitude number)




Physiological EMG activity

- No activity during relaxation
- Insertion activity
 - Caused by the mechanical stimulation by needle insertion
- Neuromuscular junction activity
 - Neuromuscular noise
 - Neuromuscular peaks






Neuromuscular noise




Neuromuscular peaks

Physiological EMG activity

- Voluntary activity
 - MU recruitment
 - Interferention
 - Willis analysis




Minimal contraction



Medium contraction 30%-40%






Maximal contraction



Pathological EMG activity


- Neuropathy
 - Slow lead, not synchronized activation causing polyphasic action potential with greater amplitude
 - Similar Mus shows relatively high frequency of firing in low, median and high effort levels
- Myopathy
 - Loss of muscle fibers,
 - Muscular dystrophy leading to dispersed action potentials with low amplitude
 - During low effort level more Mus is involved





Abnormal EMG activity




Denervated muscle fibrillation




Positive sharp peaks



Fascicle potentials




Myotonic impulses



Needle EMG Audiomyogram


Potential	Source	Sound character	Frequency	Regularity
Junction noise	MEPP	Shell	20-40Hz	Irregular
Junction potential	Terminal axon branching	Crackling	5-50Hz	Irregular
Fibrillation	Denervated muscle fiber	Rain/ ticking	0.5-10Hz	Regular
Positive sharp peak	Denervated muscle fiber	Rain/ticking	0.5-10Hz	Regular
Myotonic firing	Transversal tubulae	Motorbike start	20-150Hz	In/decreasing
Fasciculation	Motoneuron or axon	Popcorn	0.1-10Hz	Irregular
Complex repetitive firing	Denervated fibers	Engine	5-100Hz	Regular
Myokymia	Motoneuron or axon	Steps	5-60Hz	Groups
Neuromyothonia	Motoneuron or axon	Release	150-250Hz	Decreasing

Needle EMG



Stimulation EMG (Conduction study)

- Stimulation and acquisition
- Electric stimulation resulting in muscle twitch
- uncomfortable
- Examination aims:
 - Comparison of propagation speed
 - Comparison of the potential size
 - Abnormal values:
 - Nerv damage by injury, chronic pressure, inflammation, metabolic disorder and intoxication
 - Detection of lesion
 - Demyelination
 - Axonal



ADAM.



Conduction study

Muscle response on stimulus

Normal




Pathological



Direct response: M-wave

Pozdní odpovědi: F-wave

H-reflex




EMG examination

MNC	Pozice	Latence [ms]	Amplituda [mV]	CV [m/s]	Amp % [%]	Vzdálenost [mm]
Medianus-Pravá	Wrist - APB	3,9	12,0	--	--	--
Ulnaris-Pravá	Wrist - ADQ	2,5	8,0	--	--	--
Ulnaris-Levá	Wrist - ADQ	3,4	4,5	--	--	--
Medianus-Levá	Wrist - APB	4,5	6,4	--	--	--

SNC	Pozice	Latence peak [ms]	Amplituda [mV]	CV [m/s]	Vzdálenost [mm]
Median-Palm-II, Wrist-Pravá	Palm - Digit II	1,6	29,4	67	75
	Wrist - Palm	3,8	18,6	53	100
Ulnaris x Medianus Dig IV-Pravá	Ulnaris - Digit IV	3,7	9,0	67	160
	Medianus - Digit IV	3,8	10,5	52	150
Median-Palm-II, Wrist-Levá	Palm - Digit II	1,9	24,9	60	80
	Wrist - Palm	3,9	19,1	39	80
Ulnaris x Medianus Dig IV-Levá	Ulnaris - Digit IV	3,4	17,9	56	140
	Medianus - Digit IV	4,0	7,9	43	140

F-Vlna	M-Latence [ms]	M-Amplituda [mV]	Fmin [ms]	F-M [ms]	F / M [%]
Medianus-Pravá	4,0	12,0	29,0	25,1	80,0
Ulnaris-Pravá	2,4	7,8	28,4	25,9	100,0
Ulnaris-Levá	3,4	4,7	29,1	25,8	100,0
Medianus-Levá	4,5	8,2	29,4	24,9	71,0

EMG examination



EMG examination

KT Motor

Nerve / Sites	Lat ms	Amp mV	Dist cm
R MEDIANUS - APB KT			
1. Zápěstí	8,80	0,4	8
L MEDIANUS - APB KT			
1. Zápěstí	4,75	3,5	8

F Wave

Nerve	Mean F Lat ms	Min F Lat ms	Max F Lat ms	% F %	Mean F Amp mV
L TIBIAL (KNEE) - AH	52,09	51,50	52,40	100	0,2
L COMM PERONEAL - EDB	50,48	49,55	51,35	75	0,2

Tabulka EMG										
	Spontaneous					MUAP			Recruitment	Reduction
	IA	Fib	PSW	Fasc	H.F.	Amp	Dur.	PPP	Pattern	Klasifikace
L. TIB ANTERIOR	N	None	None	None		N	N	N	N	5

Závěr:

EMG vyšetření neprokázalo polyneuropatii HK či DK, ani významnější radikulopatii L5, svědčí pro středně těžkou kompresivní lézi n. medianus v oblasti zápěstí vlevo, vpravo je téměř úplná axonopatie tohoto nervu (dle anamnézy již 15 let).

MUDr. [redacted]

Pacient: [redacted]

r.č.: [redacted]

Adresa: [redacted]

Výška: 1,6 m

Age: 75 Years 9 Months

z.p.: 111

Dg.: G 61.9

Sensory NCS

Nerve / Sites	Rec. Site	Latency ms	Peak Ampl μV	Distance cm	Velocity m/s
L SURAL - Lat Malleolus					
1. Lýtko	Zevní kotník	2,00	13,5	10	50,0
R SURAL - Lat Malleolus					
1. Lýtko	Zevní kotník	2,05	23,3	10	48,8
L SUP PERONEAL - Foot					
1. Nad kotníkem	Nárt	3,10	1,5	12	38,7


KT Sens

Nerve / Sites	Rec. Site	Latency ms	Peak Ampl μV	Distance cm	Velocity m/s
L MEDIANUS - 2 prst					
1. Dlaň	II	1,85	9,8	7,3	39,5
2. Zápěstí	II	3,60	19,0	7,5	42,9
R MEDIANUS - 2 prst					
1. Dlaň	II	0,00			
2. Zápěstí	II	4,20	6,6	14,7	35,0
L MEDIANUS - x ULNAR 4, prst					
1. zápěstí - ulnar	IV	2,85	7,9	14,2	49,8
2. zápěstí - median	IV	3,50	6,0	14,2	40,6
R MEDIANUS - x ULNAR 4, prst					
1. zápěstí - ulnar	IV	2,90	9,3	14,2	49,0
2. zápěstí - median	IV	3,80	6,0	14	36,8

Motor NCS

Nerve / Sites	Lat. ms	Amp. mV	Dist. cm	RV m/s
L TIBIAL (KNEE) - AH				
1. Ankle	4,80	6,9		
2. Knee	13,30	4,5	38,2	44,9
L COMM PERONEAL - EDB				
1. Ankle	5,50	5,1		
2. Fib Head	11,90	5,0	29	45,3
3. Knee	13,55	5,1	8,5	51,5

Kineziological surface EMG




Kineziological surface EMG


Video



Muscle activity animation




Muscle activity pattern






Kineziological surface EMG

- Sequential Fourier analysis
- Window length (1 to 5 s)
- Segment overlap increases resolution
- Decision threshold for the fatigue level




Prosthetics





Other EMG applications

- Detection of myoclonic twitches
- Localisation CNS source



Other EMG applications

- Emotion detection



EMG analysis

- Artefacts
- EMG envelopes
- Feature extraction
 - Temporal
 - Frequency
- Advanced methods

EMG history



1791 - Luigi Galvani (Italy)

Relationship between electricity and muscle contractions

- **19. century**

Muscle contractions by electro stimulation

- **1907 - Louis Lapicque (France)**

Cell membrane model

- **1918 - Arthur E. Baines (England)**

First cable model of stimulus propagation

- **1928 - R. Proebster (Germany)**


First pathological signal

- **1929 – needle electrode**

následuje rychlý rozmach klinické myografie

- **2nd half of 20th century**

Electronic hand prosthesis




EMG signal


- 20-500Hz
- Minimal $F_s = 1000\text{Hz}$
- Median frequencies around 70 to 80Hz
- Do not use narrow band stopbands to suppress 50Hz hum
- It is complicated to estimate contraction level

Artefacts

ECG artefacts




Powergrid brum



50 Hz


DC offset



Nonzero DC
component


Artefacts

- Movement artefact



Electrode movement


- Amplifier saturation




clipping
 ± 0.5 V

Envelope EMG


- (band pass) EMG




- envelope EMG (threshold frequency 4 Hz)



Synchronous averaging




Mean value \pm std




Triangular detection method

Integrated EMG


Raw EMG




Integrated EMG



Integrated EMG (reset every 0.1 s)




Integrated EMG (reset at 20 uV.s)



Envelope EMG examples


Moving Average

Numerical integration



Quantitative characteristics in Frequency domain


- Median frequency
- First spectral moment
- Second spectral moment



Quantitative characteristics in Frequency domain




- Power spectral density

flexor digitorum longus (MVC) Fourier Analysis




Fatigue analysis


- erector spinae 60 seconds (50% overlap)



Zero crossing




Frequency median



First spectral moment

$$mom_1 = \frac{\sum \mathbf{f} \cdot \mathbf{I}}{\sum \mathbf{I}}$$


where \mathbf{f} is vector of corresponding frequencies.
Defines spectrum center of gravity




Second spectral moment

$$mom_2 = \sqrt{\frac{\sum \mathbf{f}^2 \cdot \mathbf{I}}{\sum \mathbf{I}} - \left(\frac{\sum \mathbf{f} \cdot \mathbf{I}}{\sum \mathbf{I}} \right)^2}$$

Is analogous to static dispersion and represent spectrum dispersion



Advanced DSP methods



Decomposition of intramuscular and surface EMG

